## ClinCal® - Calibrator

**Serum Calibrator lyophilised, Antiepileptics 4 /**

**Serum-Kalibrator lyophilisiert, Antiepileptika 4**

for HPLC Assay in Serum/Plasma

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The mean values have been assigned with the ClinRep® Complete Kit, advanced, for Antiepileptics.

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### Intended use:

ClinCal® Serum Calibrators are used for calibration of the HPLC-system. After reconstitution this lyophilised calibrator has to be prepared like a patient sample.

### Reconstitution:

Add exactly 3.0 ml HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

### Storage and stability:

This product will be stable until the expiration date when stored unopened at 2 - 8 °C. After reconstitution the stability of the analytes is:

- 12 hours at 15 - 30 °C
- 7 days at 2 - 8 °C
- 90 days at below -18 °C (avoid repeated freezing and thawing)

The variation of the filling volume (CV) is < 1 %. The average residual moisture of this lot is 0.89 %.

### Pack size:

ClinCal® Serum Calibrator Antiepileptics 4: 1 x 3 ml, **order no.: 14113**

### Caution:

The human serum which was used for manufacturing the calibrator was tested for the following infectious markers and found negative: HIV1/2- and HCV-antibodies, Hepatitis B-surface antigen, HIV- and HCV-RNA, HBV-DNA (NAT). Nevertheless the serum calibrator should be considered as potentially infectious and treated with appropriate care.

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