## ClinChek® - Control Serum Control lyophilised / Kontrollserum lyophilisiert

### ANTIEPILEPTICS / ANTIEPILEPTIKA

#### Intended use:

ClinChek® Serum Controls are used for internal quality assurance in clinical-chemical laboratories. These lyophilised controls are based on calf serum and are available with mean values in the therapeutical as well as in the elevated range. After reconstitution the controls have to be prepared like patient samples in one series of analyses.

#### Reconstitution:

Add exactly 3.0 ml of HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

#### Storage and stability:

This product will be stable until the expiration date when stored unopened at 2 - 8 °C. After reconstitution the stability of the analytes is:

- 8 hours at 15 - 30 °C (4 hours for Oxcarbazepine)
- 7 days at 2 - 8 °C
- 90 days at below -18 °C (avoid repeated freezing and thawing)

#### Notes:

The concentrations of the analytes are chosen in ranges where valid results can be obtained. According to quality assurance all ClinChek® Controls have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume (CV) is < 1 %. The average residual moisture of this lot is 1.44 %.

#### Mean values:

The mean values and confidence intervals have been established in additional independent reference laboratories with supervision of RECIPE, according to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK), with statistical methods.

#### Pack size:

ClinChek® Serum Control  
Level I, II  
2 x 5 x 3 ml, **order no.: 15082**

#### Precautions:

The calf serum which was used for manufacturing the controls should be considered as potentially infectious and therefore the controls should be treated with appropriate care.

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In Germany Phenobarbital, Phenytoin, Carbamazepine and Valproic Acid are subject to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK). For assessment of the internal quality control, please also refer to table B1a of the guideline.

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RECIPE CHEMICALS + INSTRUMENTS GmbH