## ClinChek® - Control Dried Blood Spot Control (DBS) / Trockenblut-Kontrolle (Filterspots)

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#### Intended use:

ClinChek® Dried Blood Spot Controls are used for internal quality assurance in clinical-chemical laboratories. These filter spot controls are manufactured especially for laboratories that are working in the field of neonatal screening. The controls are based on human whole blood and are available in two different ranges of concentration.

#### Pretreatment:

After having stamped a defined punching out of the blood spot, it has to be prepared like patient samples in one series of analyses. Please note that the blood volume that was used for creating one spot on the filter paper is exactly **100 µl**. The diameter of one filter spot is **18 mm**. If you use punchings of smaller diameters, you have to calculate the corresponding blood volume of it (**area ratio**).

#### Storage and stability:

Originally closed and stored below  $-18\text{ °C}$ , the dried blood spot control is stable for 36 months, but not longer than the expiration date printed on the label.

Before opening the frozen control pack, you have to temper it to room temperature, to protect against condensed water. The pack includes a pouch with desiccant (orange gel), which serves as indicator for moisture simultaneously. In the dry status the gel is coloured orange. If the capacity for water absorption is expended, the desiccant will change to blue. Then the material must be regenerated by heating for at least 5 hours at  $70\text{ °C}$  in a cabinet dryer. After cooling down, it should be enclosed to the filter cards in the aluminium bag.

Opened control packs can be stored until the expiration date, printed on the label, if they are stored dry and below  $-18\text{ °C}$ . Decreased Ornithine concentrations are possible.

#### CAUTION:

**The temperature during transport shall not exceed  $+37\text{ °C}$ . Therefore a temperature indicator is attached to the product, including appropriate instructions. Please note!**

#### Notes:

The concentrations of the analytes are chosen in relevant ranges. According to quality assurance all ClinChek® Controls have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each filter test card of one lot. The variation (CV) of the volume per blood spot is  $< 1\%$ .

#### Mean values:

The mean values and confidence intervals have been established at RECIPE, according to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK) with statistical methods. The determination of the analytes was performed by **LC-MS/MS (derivatised)**.

#### Pack size:

ClinChek® Dried Blood Spot Control (DBS)  
Level I, II  
2 x 1 x 3 spots, **order no.: MS10182**

#### Precautions:

The human whole blood which was used for manufacturing the controls was tested for the following infectious markers and found negative: HIV1/2-, HBV- and HCV-antibodies, Hepatitis B-surface antigen, HIV1- and HCV-RNA, HBV-DNA (NAT). Nevertheless, the controls should be considered as potentially infectious and treated with appropriate care.

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RECIPE CHEMICALS + INSTRUMENTS GmbH
CHEMICALS + INSTRUMENTS GmbH