# **IDK<sup>®</sup> Serotonin ELISA**

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### 14. LITERATUR

1. Badawy, A. a-B. Tryptophan: the key to boosting brain serotonin synthesis in depressive illness. *Journal of psychopharmacology (Oxford, England)* **27**, 878–93 (2013).
2. Gershon, M. D. 5-Hydroxytryptamine (serotonin) in the gastrointestinal tract. *Curr Opin Endocrinol Diabetes Obes* **20**, 14–21 (2013).
3. Gershon, M. D. Serotonin is a sword and a shield of the bowel: serotonin plays offense and defense. *Transactions of the American Clinical and Climatological Association* **123**, 268–80; discussion 280 (2012).
4. Namkung, J., Kim, H. & Park, S. Peripheral Serotonin: a New Player in Systemic Energy Homeostasis. *Molecules and Cells* **38**, 1–6 (2015).
5. Park, H.-J., Lee, S.-E., Oh, J.-H., Seo, K.-W. & Song, K.-H. Leptin, adiponectin and serotonin levels in lean and obese dogs. *BMC veterinary research* **10**, 113 (2014).
6. Reigstad, C. S., Salmonson, C. E., Rainey, J. F., Szurszewski, J. H., Linden, D. R., Sonnenburg, J. L., Farrugia, G. & Kashyap, P. C. Gut microbes promote colonic serotonin production through an effect of short-chain fatty acids on enterochromaffin cells. *The FASEB Journal* **29**, 1395–1403 (2015).

**Verwendete Symbole:**

Temperaturbegrenzung

Bestellnummer

*In-Vitro*-Diagnostikum

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Hersteller

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
# IDK<sup>®</sup> Serotonin ELISA

*For the in vitro determination of serotonin (5-HT) in human serum*

Valid from 2016-01-04

**K 6880**

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

## Table of Contents

|   |           |
|---|-----------|
| <b>1. INTENDED USE</b>                                  | <b>18</b> |
| <b>2. MATERIAL SUPPLIED</b>                             | <b>18</b> |
| <b>3. MATERIAL REQUIRED BUT NOT SUPPLIED</b>            | <b>18</b> |
| <b>4. PREPARATION AND STORAGE OF REAGENTS</b>           | <b>19</b> |
| <b>5. PREPARATION AND STORAGE OF SAMPLES</b>            | <b>20</b> |
| <b>6. ASSAY PROCEDURE</b>                               | <b>20</b> |
| <i>Principle of the test</i>                            | 20        |
| <i>Sample preparation procedure</i>                     | 21        |
| <i>Test procedure</i>                                   | 22        |
| <b>7. RESULTS</b>                                       | <b>23</b> |
| <b>8. LIMITATIONS</b>                                   | <b>24</b> |
| <b>9. QUALITY CONTROL</b>                               | <b>25</b> |
| <i>Reference range</i>                                  | 25        |
| <b>10. PERFORMANCE CHARACTERISTICS</b>                  | <b>25</b> |
| <i>Precision and reproducibility</i>                    | 25        |
| <i>Spiking recovery</i>                                 | 26        |
| <i>Dilution recovery</i>                                | 26        |
| <i>Analytical sensitivity</i>                           | 26        |
| <i>Specificity</i>                                      | 27        |
| <b>11. PRECAUTIONS</b>                                  | <b>27</b> |
| <b>12. TECHNICAL HINTS</b>                              | <b>27</b> |
| <b>13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b> | <b>28</b> |
| <b>14. REFERENCES</b>                                   | <b>28</b> |

## 1. INTENDED USE

This Immundiagnostik assay is intended for the quantitative determination of serotonin (5-hydroxytryptamine, 5-HT) in human serum. For *in vitro* diagnostic use only.

## 2. MATERIAL SUPPLIED

| Catalog No.      | Label            | Kit Components  | Quantity       |
|------------------|------------------|---|----------------|
| K 6880           | PLATE            | Holder with precoated strips                                | 12 x 8 wells   |
| K 6880           | STD              | Standards, ready to use<br>(0, 15, 50, 100, 250, 750 ng/ml) | 6 x 200 µl     |
| K 6880<br>K 6880 | CTRL 1<br>CTRL 2 | Controls, ready to use<br>(see specification for range)     | 2 x 200 µl     |
| K 6880           | WASHBUF          | ELISA wash buffer concentrate, 10x                          | 2 x 100 ml     |
| K 6880           | CONJ             | Conjugate (peroxidase-labeled),<br>concentrate              | 1 x 50 µl      |
| K 6880           | CONJBUF          | Conjugate stabilizing buffer, ready to use                  | 1 x 13 ml      |
| K 6880           | REABUF           | Reaction buffer, ready to use                               | 2 x 30 ml      |
| K 6880           | DER              | Derivatization reagent                                      | 2 x 50 mg      |
| K 6880           | DMF              | Dimethylformamide (DMF)                                     | 1 x 7 ml       |
| K 6880           | ASYREAG          | Assay reagent, lyophilized                                  | 2 x 30 ml lyo. |
| K 6880           | SUB              | TMB substrate (tetramethylbenzidine),<br>ready to use       | 1 x 15 ml      |
| K 6880           | STOP             | ELISA stop solution, ready to use                           | 1 x 15 ml      |

For reorders of single components, use the catalogue number followed by the label as product number.

## 3. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water\*
- Calibrated precision pipets and 10-1000 µl tips

- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Vortex
- Centrifuge, 3000 x *g*
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 6)

\* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥18.2 MΩ cm).

#### 4. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than once ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 2 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- Dilute the **wash buffer concentrate (WASHBUF)** with ultra pure water **1:10** before use (**100 ml WASHBUF + 900 ml ultra pure water**), mix well. Crystals may occur due to high salt concentration in the stock solution. The crystals must be redissolved at room temperature or at 37 °C using a water bath before dilution. The buffer concentrate is stable at **2-8 °C** until the expiry date stated on the label. Diluted buffer solution can be stored in a closed flask at **2-8 °C for one month.**
- Store **Standards (STD)** and **controls (CTRL1, CTRL2)** frozen at **-20 °C**, thaw before use in the test and mix thoroughly for several seconds on a vortex mixer. Any occurring solid matter can be dissolved by gentle warming. Re-freeze after use. Standards and controls can be thawed and refrozen for several times.
- Store **reaction buffer (REABUF)** frozen at **-20 °C**. Bring to room temperature and mix before use in the test. Re-freeze after use.
- Dissolve the content of one vial of **derivatization reagent (DER) (50 mg) in 3 ml DMF**. CAUTION! DMF is toxic (see chapter 11 - precautions). Put the vial on a horizontal shaker for 5 min. DER must be **prepared immediately**

**before use.** When more than one vial is to be used, combine the contents and mix prior to use. Dispose of any rest of the reagent after use. The ELISA kit can be separated into two performances by providing two DER vials. Please note: DMF attacks all plastics but not polypropylene products and laboratory glass.

- Dissolve the **assay reagent (ASYREAG)** in **30 ml of diluted wash buffer (WASHBUF)**, mix thoroughly. Put the vial on a horizontal shaker for 10 min and ensure complete dissolving of any solid matter. When more than one vial is to be used, combine the contents in a separate container and mix prior to use. Discard any rest after use. The ELISA kit can be separated into two performances by providing two vials of assay buffer
- Dilute the **peroxidase conjugate (CONJ) 1:501** with conjugate stabilizing buffer (CONJBUF) (e.g. **25 µl CONJ + 12 ml CONJBUF**, prepare only the required amount). The undiluted POD conjugate is stable at **2-8 °C** until the expiry date stated on the label. Diluted POD conjugate can be stored at **2-8 °C for 1 week.**
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2-8 °C.**

## 5. PREPARATION AND STORAGE OF SAMPLES

- The serum samples are stable for up to 48 hours at 2-8 °C. For longer storage keep samples frozen at -20 °C. As serotonin is sensitive to light, store samples in a cool and dark place after collection.
- Lipemic or hemolytic samples may give erroneous results and should not be used for analysis.
- Samples with visible amounts of precipitates should be centrifuged.
- The serum samples are analyzed **undiluted.**
- For sample preparation a derivatization reagent (DER) for derivatization of serotonin is added (details are given in the sample preparation procedure).

## 6. ASSAY PROCEDURE

### *Principle of the test*

This assay is based on the method of competitive enzyme linked immunoassays. The sample preparation includes the addition of a

derivatization reagent for serotonin derivatization. Assay reagent is also added containing serotonin-derivative (tracer). Afterwards, the treated samples are incubated in wells of a microtiter plate coated with a polyclonal antibody against serotonin-derivative. During the incubation period, the target serotonin in the sample competes with the tracer for the binding of the polyclonal antibodies immobilized on the wall of the microtiter wells. The serotonin in the sample displaces the tracer out of the binding to the antibodies. Therefore, the concentration of the antibody-bound tracer is inverse proportional to the serotonin concentration in the sample.

During the second incubation step, a peroxidase conjugate is added to each microtiter well to detect the tracer. After washing away the unbound components tetramethylbenzidine (TMB) is added as a peroxidase substrate. Finally, the enzymatic reaction is terminated by an acidic stop solution. The color changes from blue to yellow, and the absorbance is measured in a photometer at 450 nm. The intensity of the yellow color is inverse proportional to the serotonin concentration in the sample; this means, high serotonin concentration in the sample reduces the concentration of antibody-bound tracer and lowers the photometric signal.

A dose response curve of absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standards. serotonin present in the patient samples is determined directly from this curve.

### *Sample preparation procedure*

Derivatization of standards (STD), controls (CTRL) and samples (SAMPLE) is carried out in single analysis in vials (e.g. **2 ml** vials).

The reagents provided with this kit are sufficient for up to 48 derivatizations, which are transferred in duplicate determinations to the wells of the microtiter plate.

1. Bring all reagents and samples to **room temperature** (15-30 °C) and mix well.
2. Add **25 µl of standards (STD)**, **25 µl of controls (CTRL)**, and **25 µl of samples (SAMPLE)** in the corresponding vials.
3. Add **500 µl of reaction buffer (REABUF)** into each vial (STD, CTRL, SAMPLE).

4. Add **50 µl** of freshly prepared **derivatization reagent (DER)** into each vial (STD, CTRL, SAMPLE), **mix thoroughly** by repeated inversion or several seconds on a vortex mixer and incubate for **30 min at room temperature** (15-30°C) on a horizontal **shaker** (180-240 rpm).
5. Afterwards add **1 ml of assay reagent (ASYREAG)** into each vial and mix well.

**2 x 100 µl of each treated sample (STD, CTRL, SAMPLE)** are used in the ELISA as duplicates.

### *Test procedure*

6. Mark the positions of standards (STD)/ controls (CTRL)/samples (SAMPLE) in duplicate on a **protocol sheet**.
7. Take as many **microtiter plate strips (PLATE)** as needed from kit. Store unused strips covered at 2-8°C. Strips are stable until the expiry date stated on the label.
8. Wash each well **5 x** with **250 µl of diluted wash buffer**. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
9. For the analysis in duplicate take **2 x 100 µl** of the **derivatized standards/ controls/ samples (STD/ CTRL/ SAMPLE)** out of the vials and add into the respective wells.
10. Cover the plate and incubate for **1 hour at room temperature** (15-30 °C) on a horizontal shaker (180-240 rpm).
11. Discard the contents of each well. Wash each well **5 x** with **250 µl of diluted wash buffer**. After the final washing step the inverted microtiter plate should be firmly tapped on absorbent paper.
12. Add **100 µl** of diluted **peroxidase conjugate (CONJ)** into each well.
13. Cover the plate and incubate for **30 min at room temperature** (15-30 °C) on a horizontal **shaker** (180-240 rpm).

14. Discard the contents of each well. Wash each well **5 x** with **250 µl** of **diluted wash buffer**. After the final washing step the inverted microtiter plate should be firmly tapped on absorbent paper.
15. Add **100 µl of TMB substrate (SUB)** into each well.
16. Incubate for **7-12 min** at **room temperature** (15-30 °C) in the dark\*.
17. Add **100 µl of stop solution (STOP)** into each well, mix thoroughly.
18. Determine **absorption immediately** with an ELISA reader at **450 nm** against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at **405 nm** against 620 nm (690 nm) as a reference.

\* The intensity of the color change is temperature sensitive. We recommend to observe the color change and to stop the reaction upon good differentiation.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

## 7. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the "4 parameter algorithm".

### 1. 4 parameter algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e.g. 0.001).

### 2. Point-to-point calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

### 3. Spline algorithm

We recommend a linear ordinate for optical density and a linear abscissa for concentration.


The plausibility of the duplicate values should be examined before automatically evaluating the results. If this option is not available within the used program, the duplicate values should be evaluated manually.

#### Serum

The concentrations can be determined directly from the standard curve in ng/ml. **No factor** is needed.

If a sample has been diluted, multiply the obtained result by the dilution factor used.

In the following, an example of a standard curve is given; do not use it for the calculation of your results.

## 8. LIMITATIONS

Samples with an OD lower than the OD of the highest standard should be diluted with reaction buffer (REABUF) and re-assayed. Please consider this dilution factor when calculating the results.

## 9. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analyzed with each run. Results generated from the analysis of control samples should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control samples are outside of the acceptable limits.

### *Reference range*

Based on internal studies with serum samples of apparently healthy persons (n = 36) a mean value of 225 ng/ml was calculated. The standard deviation was 86 ng/ml.

**Serum mean value  $\pm$  2 x standard deviation:** **225  $\pm$  172 ng/ml**

**Normal range:** **53 – 397 ng/ml**

We recommend each laboratory to establish its own reference range.

## 10. PERFORMANCE CHARACTERISTICS

### *Precision and reproducibility*

#### **Intra-assay (n = 8)**

| Sample | Serotonin [ng/ml] | CV [%] |
|--------|-------------------|--------|
| 1      | 171.5             | 8.4    |
| 2      | 80.3              | 12.1   |

#### **Inter-assay (n = 6)**

| Sample | Serotonin [ng/ml] | CV [%] |
|--------|-------------------|--------|
| 1      | 293.1             | 9.0    |
| 2      | 400.9             | 18.9   |

### Spiking recovery

Two serum samples were spiked with different serotonin concentrations and measured in this assay. The mean recovery rate was 102.5 % (n = 2).

| Sample [ng/ml] | Spike [ng/ml] | Serotonin expected [ng/ml] | Serotonin measured [ng/ml] | Recovery [%] |
|----------------|---------------|----------------------------|----------------------------|--------------|
| 180.6          | 100           | 280.6                      | 298.4                      | 106.3        |
| 180.6          | 200           | 380.6                      | 402.3                      | 105.7        |
| 100.3          | 100           | 200.3                      | 186.3                      | 93.0         |
| 100.3          | 200           | 300.3                      | 315.6                      | 105.1        |

### Dilution recovery

Two spiked serum samples were diluted with reaction buffer. The mean recovery was 99.4 % (n = 2).

| Sample [ng/ml] | Dilution | Serotonin expected [ng/ml] | Serotonin measured [ng/ml] | Recovery [%] |
|----------------|----------|----------------------------|----------------------------|--------------|
| 175.1          | 1:2      | 87.6                       | 92.5                       | 105.6        |
| 175.1          | 1:4      | 43.8                       | 48.0                       | 109.5        |
| 246.8          | 1:2      | 123.4                      | 114.0                      | 92.3         |
| 246.8          | 1:4      | 61.7                       | 55.7                       | 90.2         |

### Analytical sensitivity

The zero-standard was measured 36 times. The detection limit was set as  $B_0 - 2 \text{ SD}$  and estimated to be 2.3 ng/ml.

| Sample        | mean value [OD] | 2 x standard deviation (2 x SD) | Detection limit [ng/ml] |
|---------------|-----------------|---------------------------------|-------------------------|
| Zero-standard | 2.186           | 0.13                            | 2.3                     |

## Specificity

Specificity was tested by measuring the cross-reactivity against compounds with structural similarity to serotonin. The specificity is calculated in percent in relation to the serotonin binding activity.

|                       |          |
|-----------------------|----------|
| 3-Indoleacrylic acid  | < 0.01%  |
| Indole-3-pyruvic acid | < 0.04%  |
| 3-Indoleacetic acid   | < 0.05 % |
| 5-Methoxytryptophol   | < 0.08 % |
| L-5-OH-tryptophan     | < 0.01 % |

## 11. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulfuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breathe vapour and avoid inhalation.
- DMF is toxic. It may cause harm to the unborn child and is harmful by inhalation and in contact with skin. Work under hood and apply preventive skin protection. In case of skin or eye contact flush with plenty of water and get medical attention.

## 12. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore, we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.

- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

### 13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- *IDK®* is a trademark of Immundiagnostik AG.
- Incubation time, incubation temperature, and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

### 14. REFERENCES

1. Badawy, A. a-B. Tryptophan: the key to boosting brain serotonin synthesis in depressive illness. *Journal of psychopharmacology (Oxford, England)* **27**, 878–93 (2013).
2. Gershon, M. D. 5-Hydroxytryptamine (serotonin) in the gastrointestinal tract. *Curr Opin Endocrinol Diabetes Obes* **20**, 14–21 (2013).
3. Gershon, M. D. Serotonin is a sword and a shield of the bowel: serotonin plays offense and defense. *Transactions of the American Clinical and Climatological Association* **123**, 268–80; discussion 280 (2012).
4. Namkung, J., Kim, H. & Park, S. Peripheral Serotonin: a New Player in Systemic Energy Homeostasis. *Molecules and Cells* **38**, 1–6 (2015).

5. Park, H.-J., Lee, S.-E., Oh, J.-H., Seo, K.-W. & Song, K.-H. Leptin, adiponectin and serotonin levels in lean and obese dogs. *BMC veterinary research* **10**, 113 (2014).
6. Reigstad, C. S., Salmonson, C. E., Rainey, J. F., Szurszewski, J. H., Linden, D. R., Sonnenburg, J. L., Farrugia, G. & Kashyap, P. C. Gut microbes promote colonic serotonin production through an effect of short-chain fatty acids on enterochromaffin cells. *The FASEB Journal* **29**, 1395–1403 (2015).

**Used symbols:**

|   |                                    |   |                                   |
|---|------------------------------------|---|-----------------------------------|
|    | Temperature limitation             |    | Catalogue Number                  |
|   | In Vitro Diagnostic Medical Device |   | To be used with                   |
|  | Manufacturer                       |  | Contains sufficient for <n> tests |
|  | Lot number                         |  | Use by                            |