IDK[®] Tryptophan ELISA

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15. LITERATUR

1. Brandacher, G. et al., 2007. Chronic immune activation underlies morbid obesity: is IDO a key player? *Current drug metabolism*, **8**(3), pp.289–95.
2. Miura, H. et al., 2008. A link between stress and depression: shifts in the balance between the kynurenine and serotonin pathways of tryptophan metabolism and the etiology and pathophysiology of depression. *Stress (Amsterdam, Netherlands)*, **11**(3), pp.198–209.
3. Wilson, S.T. et al., 2009. The tryptophan hydroxylase-1 A218C polymorphism is associated with diagnosis, but not suicidal behavior, in borderline personality disorder. *American journal of medical genetics. Part B, Neuropsychiatric genetics : the official publication of the International Society of Psychiatric Genetics*, **150B**(2), pp.202–8.
4. Suzuki, Y. et al., 2012. Serum indoleamine 2,3-dioxygenase activity predicts prognosis of pulmonary tuberculosis. *Clinical and vaccine immunology : CVI*, **19**(3), pp.436–42.
5. Kim, H. et al., 2012. Brain indoleamine 2,3-dioxygenase contributes to the comorbidity of pain and depression. *The Journal of clinical investigation*, **122**(8), pp.2940–54.
6. Gupta, N.K. et al., 2012. Serum analysis of tryptophan catabolism pathway: correlation with Crohn's disease activity. *Inflammatory bowel diseases*, **18**(7), pp.1214–20.
7. Sulo, G. et al., 2013. Neopterin and kynurenine-tryptophan ratio as predictors of coronary events in older adults, the Hordaland Health Study. *International journal of cardiology*, **168**(2), pp.1435–40.

8. Pedersen, E.R. et al., 2013. Urinary excretion of kynurenine and tryptophan, cardiovascular events, and mortality after elective coronary angiography. *European heart journal*, **34**(34), pp.2689–96.
9. Creelan, B.C. et al., 2013. Indoleamine 2,3-dioxygenase activity and clinical outcome following induction chemotherapy and concurrent chemoradiation in Stage III non-small cell lung cancer. *Oncoimmunology*, **2**(3), p.e23428.
10. Ciorba, M.A., 2013. Indoleamine 2,3 dioxygenase in intestinal disease. *Current opinion in gastroenterology*, **29**(2), pp.146–52.
11. Ristagno, G. et al., 2014. Early activation of the kynurenine pathway predicts early death and long-term outcome in patients resuscitated from out-of-hospital cardiac arrest. *Journal of the American Heart Association*, **3**(4), pp.e001094–e001094.
12. Grozdics, E. et al., 2014. B7 costimulation and intracellular indoleamine-2,3-dioxygenase (IDO) expression in peripheral blood of healthy pregnant and non-pregnant women. *BMC pregnancy and childbirth*, **14**(1), p.306.
13. Dolina, S. et al., 2014. Attention-deficit hyperactivity disorder (ADHD) as a pyridoxine-dependent condition: urinary diagnostic biomarkers. *Medical hypotheses*, **82**(1), pp.111–6.
14. Chuang, S.-C. et al., 2014. Circulating biomarkers of tryptophan and the kynurenine pathway and lung cancer risk. *Cancer epidemiology, biomarkers & prevention : a publication of the American Association for Cancer Research, cosponsored by the American Society of Preventive Oncology*, **23**(3), pp.461–8.
15. Yan, E.B. et al., 2015. Activation of the kynurenine pathway and increased production of the excitotoxin quinolinic acid following traumatic brain injury in humans. *Journal of neuroinflammation*, **12**(1), p.110.

Verwendete Symbole:

Temperaturbegrenzung

Bestellnummer

*In-Vitro*-Diagnostikum

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Hersteller

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
Achtung

IDK[®] Tryptophan ELISA

*For the in vitro determination of L-tryptophan in
human dried blood samples*

Valid from 2017-01-04

REF K 7730DBS

IVD **CE**

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com www.immundiagnostik.com

Table of Contents

1. INTENDED USE	17
2. INTRODUCTION	17
3. MATERIAL SUPPLIED	18
4. MATERIAL REQUIRED BUT NOT SUPPLIED	18
5. STORAGE AND PREPARATION OF REAGENTS	19
6. PREPARATION OF SAMPLES	20
7. ASSAY PROCEDURE	20
<i>Principle of the test</i>	20
<i>Sample preparation procedure</i>	21
<i>Test procedure</i>	21
8. RESULTS	23
9. LIMITATIONS	23
10. QUALITY CONTROL	23
<i>Reference range</i>	24
11. PERFORMANCE CHARACTERISTICS	24
<i>Analytical Sensitivity</i>	24
<i>Precision and reproducibility</i>	24
12. PRECAUTIONS	24
13. TECHNICAL HINTS	25
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	25
15. REFERENCES	26

1. INTENDED USE

This Immundiagnostik assay is intended for the quantitative determination of L-tryptophan in human dried blood samples. For *in vitro* diagnostic use only.

2. INTRODUCTION

Tryptophan is an essential amino acid which cannot be synthesized by humans, and therefore must be supplied in diet. Tryptophan plays an important role in the synthesis of serotonin, as it is the only source for synthesizing serotonin (see figure 1).


Figure 1: Synthesis of serotonin and melatonin via tryptophan (left), and shifting of the metabolism of tryptophan to the kynurenine pathway via the enzyme indoleamine 2,3-dioxygenase (right).

Serotonin is released from the appendices of a widespread transmitter system with global modulating effects on different areas of life. The activity of this system depends not least on the actual tryptophan level.

When tryptophan in plasma decreases, it is not available for serotonin synthesis and the following melatonin production. As a result, fluctuations of mood, aggression, sleep and eating disorder, or depression can occur. The symptoms of these diseases can be eased by administration of tryptophan. However, such a treatment requires the knowledge of the tryptophan level, which can be determined with this test.

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 7730DBS	PLATE	Holder with precoated strips	12 x 8 wells
K 7730DBS	STD	Standards, lyophilised (see specification for concentrations)	1 x 6 vials
K 7730DBS	CTRL 1	Control, lyophilised (see specification for range)	1 x 1 vial
K 7730DBS	CTRL 2	Control, lyophilised (see specification for range)	1 x 1 vial
K 7730DBS	WASHBUF	ELISA wash buffer concentrate, 10x	2 x 100 ml
K 7730DBS	AB	L-tryptophan antibody, lyophilised	4 x 1 vial
K 7730DBS	CONJ	Conjugate concentrate, peroxidase-labelled	1 x 65 µl
K 7730DBS	CONJBUF	Conjugate stabilizing buffer, ready-to-use	1 x 13 ml
K 7730DBS	REABUF	Reaction buffer, ready-to-use	2 x 28 ml
K 7730DBS	DER	Lyophilised derivatisation reagent	4 x 25 mg
K 7730DBS	DMSO	Dimethylsulfoxide (DMSO)	1 x 7 ml
K 7730DBS	SUB	TMB substrate (tetramethylbenzidine), ready-to-use	1 x 15 ml
K 7730DBS	STOP	ELISA stop solution, ready-to-use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Calibrated precision pipettors and 10–1000 µl tips
- Foil to cover the microtiter plate

- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Vortex
- Centrifuge, 3000 g
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** has to be diluted with ultra pure water **1:10** before use (100 ml WASHBUF + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or in a water bath at 37 °C before dilution of the buffer solutions. The **WASHBUF** is stable at **2–8 °C** until the expiry date stated on the label. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2–8 °C for one month**.
- The **lyophilised standards (STD) and controls (CTRL)** are stable at **2–8 °C** until the expiry date stated on the label. Before use, reconstitute the standards and controls with **200 µl reaction buffer (REABUF)**. Mix thoroughly and allow to dissolve for 15 min on a horizontal shaker. Standards and controls (reconstituted STD and CTRL) can be stored at **-20 °C**.
- **DMSO** crystallises at 4 °C. Before use, dissolve the crystals at room temperature or in a water bath.
- Dissolve the content of one vial of **lyophilised derivatization reagent (DER) (25 mg)** in **1.5 ml DMSO**. Put the vial on a horizontal shaker for 5 min. **Derivatization reagent** (reconstituted DER) must be **prepared immediately before use**. When more than one vial is to be used, combine the contents and mix prior to use. Discard any rest of the reagent after use. **Please note:** DMSO attacks all plastics but not polypropylene products and laboratory glass.

- Reconstitute the content of **one vial of lyophilised L-tryptophan antibody (AB)** in **3 ml wash buffer**. When more than one vial is to be used, combine the contents and mix prior to use. **L-tryptophan antibody** (reconstituted AB) can be stored at **-20 °C for 1 month**.
- Dilute the **peroxidase conjugate concentrate (CONJ) 1:201 with conjugate stabilizing buffer (CONJBUF)** (e.g. 60 µl CONJ + 12 ml CONJBUF, prepare only the required amount). The CONJ is stable at 2–8 °C until the expiry date stated on the label. **Conjugate** (1:201 diluted CONJ) can be stored at **2–8 °C for 1 week**.
- All other test reagents are ready-to-use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8 °C**.

6. PREPARATION OF SAMPLES

Samples are used undiluted.

We recommend DrySpot-ID (catalogue no DZ9020ID or DZ9021ID) as dried blood spot carrier. The moistened cards are stable for 6 days at room temperature.

Place the filter paper containing the dried sample in a 1.5 ml tube. Add **500 µl reaction buffer** (REABUF) and mix well. After incubation for 30 minutes at room temperature (15–30 °C) mix again well.

For sample preparation, a derivatisation reagent (DER) for derivatisation of L-tryptophan is added (details are given in the sample preparation procedure).

7. ASSAY PROCEDURE

Principle of the test

This assay is based on the method of competitive enzyme linked immunoassays. The sample preparation includes the addition of a derivatization reagent for L-tryptophan derivatization. Afterwards, the treated samples and a polyclonal L-tryptophan antiserum are incubated in the wells of a microtiter plate coated with L-tryptophan derivative (tracer).

During the incubation period, the target L-tryptophan in the sample competes with the tracer immobilized on the wall of the microtiter wells for the binding of the polyclonal antibodies. The L-tryptophan in the sample displaces the antibodies out of the binding to the tracer. Therefore the concentration of the tracer-bound antibody is inverse proportional to the L-tryptophan concentration in the sample.

During the second incubation step a peroxidase-conjugated antibody is added to each microtiter well to detect the L-tryptophan antibodies. After washing away the unbound components tetramethylbenzidine (TMB) is added as a peroxidase substrate. Finally, the enzymatic reaction is terminated by an acidic stop solution. The color changes from blue to yellow and the absorbance is measured in a photometer at 450 nm. The intensity of the yellow color is inverse proportional to the L-tryptophan concentration in the sample; this means, high L-tryptophan concentration in the sample reduces the concentration of tracer-bound antibodies and lowers the photometric signal. A dose response curve of absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standard. L-tryptophan present in the patient samples is determined directly from this curve.

Sample preparation procedure

Derivatisation of standards, controls and samples is carried out in single analysis in vials (e.g. 1.5 ml vials).

We recommend to prepare one derivatisation per standard, sample and control and pipet it in duplicate into the microtitre plate.

1.	Bring all reagents and samples to room temperature (15–30°C) and mix well.
2.	Add each 25 µl standards or controls into the corresponding vials.
3.	Add 500 µl reaction buffer (REABUF) to the vials of standards and controls.
4.	Add 50 µl of freshly prepared derivatization reagent (reconstituted DER) into each vial (STD, CTRL, dried blood samples) and mix thoroughly, e.g. by repeated inversion or several seconds on a vortex mixer. Incubate for 45 min at room temperature (15–30°C) on a horizontal shaker .
5.	Centrifuge the vials containing the dried blood samples for 5 min at 3000 g. Use the supernatant for the following procedures.

2 x 25 µl of each treated sample (standards, controls, samples) are used in the ELISA as duplicates.

Test procedure

Bring all **reagents and samples to room temperature** (15–30°C) and mix well.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

6.	Mark the positions of standards/controls/samples on a protocol sheet .
7.	Take as many microtiter strips as needed from kit. Store unused strips covered at 2–8 °C. Strips are stable until expiry date stated on the label. Please note: do not wash the microtiter strips!
8.	Add 25 µl of the derivatised standards/controls/ samples into the respective wells.
9.	Add 100 µl L-tryptophan antibody into each well.
10.	Cover the wells and incubate for 2 hours at room temperature (15–30 °C) on a horizontal shaker .
11.	Discard the contents of each well and wash 5 times with 250 µl wash buffer . After the final washing step, remove residual wash buffer by firmly tapping the microtiter strips on absorbent paper..
12.	Add 100 µl conjugate into each well.
13.	Cover the wells and incubate for 1 hour at room temperature (15–30 °C) on a horizontal shaker .
14.	Discard the contents of each well and wash 5 times with 250 µl wash buffer . After the final washing step, remove residual wash buffer by firmly tapping the microtiter strips on absorbent paper..
15.	Add 100 µl TMB substrate into each well.
16.	Incubate for 10–15 min at room temperature (15–30 °C) in the dark* .
17.	Add 100 µl ELISA stop solution (STOP) into each well, mix thoroughly.
18.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (595 nm, 630nm, 650 nm or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm (595 nm, 630nm, 650 nm or 690 nm) as a reference.

* The intensity of the color change is temperature sensitive. We recommend observing the color change and stopping the reaction upon good differentiation.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the "4 parameter algorithm".

1. 4 parameter algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e.g. 0.001).

2. Point-to-point calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

3. Spline algorithm

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the programme used, the duplicate values should be evaluated manually.

Dried blood samples

The obtained results have to be multiplied with the **factor 2** to get the actual concentrations.

In case **another dilution factor** has been used, multiply the obtained result with the dilution factor used.

9. LIMITATIONS

Samples with an OD higher than the OD of the highest standard should be further diluted and re-assayed. For the following analysis, the changed dilution factor has to be taken into consideration.

10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the

same assay one or more values of the quality control sample are outside the acceptable limits.

Reference range

Based on internal studies with serum samples of apparently healthy persons a mean value of 63.5 µmol/l was estimated (n = 70). The standard deviation was 10.1 µmol/l.

Mean value ± 2x standard deviation: 63.5 ± 20.2 µmol/l

Normal range: 43.3–83.7 µmol/l

We recommend each laboratory to establish its own reference range.

11. PERFORMANCE CHARACTERISTICS

Analytical Sensitivity

Limit of blank, LoB

2.7 µmol/l

The evaluation was performed according to the CLSI guideline EP-17-A2. The specified accuracy goal for the LoQ was 20% CV.

Precision and reproducibility

Intra-Assay (n = 20)

Sample	L-tryptophan [µmol/l]	CV [%]
1	42.30	4.95
2	60.30	8.98

Inter-Assay (n = 20)

Sample	L-tryptophan [µmol/l]	CV [%]
1	75.3	10.3
2	106.39	10.1

12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.

- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.


- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immunodiagnostik AG along with a written complaint.

15. REFERENCES

1. Brandacher, G. et al., 2007. Chronic immune activation underlies morbid obesity: is IDO a key player? *Current drug metabolism*, **8**(3), pp.289–95.
2. Miura, H. et al., 2008. A link between stress and depression: shifts in the balance between the kynurenine and serotonin pathways of tryptophan metabolism and the etiology and pathophysiology of depression. *Stress (Amsterdam, Netherlands)*, **11**(3), pp.198–209.
3. Wilson, S.T. et al., 2009. The tryptophan hydroxylase-1 A218C polymorphism is associated with diagnosis, but not suicidal behavior, in borderline personality disorder. *American journal of medical genetics. Part B, Neuropsychiatric genetics : the official publication of the International Society of Psychiatric Genetics*, **150B**(2), pp.202–8.
4. Suzuki, Y. et al., 2012. Serum indoleamine 2,3-dioxygenase activity predicts prognosis of pulmonary tuberculosis. *Clinical and vaccine immunology : CVI*, **19**(3), pp.436–42.
5. Kim, H. et al., 2012. Brain indoleamine 2,3-dioxygenase contributes to the comorbidity of pain and depression. *The Journal of clinical investigation*, **122**(8), pp.2940–54.
6. Gupta, N.K. et al., 2012. Serum analysis of tryptophan catabolism pathway: correlation with Crohn's disease activity. *Inflammatory bowel diseases*, **18**(7), pp.1214–20.
7. Sulo, G. et al., 2013. Neopterin and kynurenine-tryptophan ratio as predictors of coronary events in older adults, the Hordaland Health Study. *International journal of cardiology*, **168**(2), pp.1435–40.
8. Pedersen, E.R. et al., 2013. Urinary excretion of kynurenine and tryptophan, cardiovascular events, and mortality after elective coronary angiography. *European heart journal*, **34**(34), pp.2689–96.
9. Creelan, B.C. et al., 2013. Indoleamine 2,3-dioxygenase activity and clinical outcome following induction chemotherapy and concurrent chemoradiation in Stage III non-small cell lung cancer. *Oncoimmunology*, **2**(3), p.e23428.
10. Ciorba, M.A., 2013. Indoleamine 2,3 dioxygenase in intestinal disease. *Current opinion in gastroenterology*, **29**(2), pp.146–52.

11. Ristagno, G. et al., 2014. Early activation of the kynurenine pathway predicts early death and long-term outcome in patients resuscitated from out-of-hospital cardiac arrest. *Journal of the American Heart Association*, **3**(4), pp.e001094–e001094.
12. Grozdics, E. et al., 2014. B7 costimulation and intracellular indoleamine-2,3-dioxygenase (IDO) expression in peripheral blood of healthy pregnant and non-pregnant women. *BMC pregnancy and childbirth*, **14**(1), p.306.
13. Dolina, S. et al., 2014. Attention-deficit hyperactivity disorder (ADHD) as a pyridoxine-dependent condition: urinary diagnostic biomarkers. *Medical hypotheses*, **82**(1), pp.111–6.
14. Chuang, S.-C. et al., 2014. Circulating biomarkers of tryptophan and the kynurenine pathway and lung cancer risk. *Cancer epidemiology, biomarkers & prevention : a publication of the American Association for Cancer Research, cosponsored by the American Society of Preventive Oncology*, **23**(3), pp.461–8.
15. Yan, E.B. et al., 2015. Activation of the kynurenine pathway and increased production of the excitotoxin quinolinic acid following traumatic brain injury in humans. *Journal of neuroinflammation*, **12**(1), p.110.

Used symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		