ClinCal[®] - Calibrator

Plasma Calibrator lyophilised, for Homocysteine /

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The mean value has been assigned with the ClinRep[®] Complete Kit for Homocysteine. The uncertainty of the mean value is $\pm 1.4\%$ at a confidence level of 95 %. The mean value is traceable to **NIST-SRM1955**.

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Intended use:

ClinCal[®] Plasma Calibrators are used for calibration of the HPLC-system. After reconstitution this lyophilised calibrator has to be prepared like a patient sample.

Reconstitution:

Add exactly 3.0 ml HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

Storage and stability:

Originally closed and stored at 2 - 8 °C, the lyophilised plasma calibrator is stable for 48 months, but not longer than the expiration date printed on the label. After reconstitution the stability of the analytes is:

- at least 12 hours when stored at 15 - 30 °C
- at least 7 days when stored at 2 - 8 °C
- at least 3 months when stored below -18 °C
(avoid repeated freezing and thawing)

The variation of the filling volume (CV) is $< 1\%$. The average residual moisture of this lot is 1.14 %.

Pack size:

ClinCal[®] Plasma Calibrator for Homocysteine: 1 x 3 ml, **order no.: 23013**

Caution:

The human plasma which was used for manufacturing the calibrator was tested for the following infectious markers and found negative: HIV1/2-, HBV- and HCV-antibodies, Hepatitis B-surface antigen, HIV- and HCV-RNA, HBV-DNA (NAT). Nevertheless, the plasma calibrator should be considered as potentially infectious and treated with appropriate care.

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