ClinCal® - Calibrator

Plasma Calibrator, lyophilised, for β -Carotene /

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Intended use:

ClinCal® Plasma Calibrators are used for calibration of the HPLC-system. After reconstitution this lyophilised calibrator has to be prepared like a patient sample.

Reconstitution:

Add exactly 0.5 ml HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

Storage and stability:

Originally closed and stored below -18 °C in the dark, the lyophilised plasma calibrator is stable for 24 months, but not longer than the expiration date printed on the label. After reconstitution the stability of the analyte is:

- at least 1 day when stored at 15 - 30 °C in the dark
- at least 7 days when stored at 2 - 8 °C in the dark
- at least 3 months when stored below -18 °C in the dark
(avoid repeated freezing and thawing)

The variation of the filling volume (CV) is < 1 %. The average residual moisture of this lot is 0.87 %.

Pack size:

ClinCal® Plasma Calibrator for β -Carotene: 5 x 0.5 ml, **order no.: 30013**

Caution:

The human plasma which was used for manufacturing the calibrator was tested for the following infectious markers and found negative: HIV1/2-, HBV- and HCV-antibodies, Hepatitis B-surface antigen, HIV- und HCV-RNA, HBV-DNA (NAT). Nevertheless, the plasma calibrator should be considered as potentially infectious and treated with appropriate care.

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