ClinTest® - Standard

Standard Solution for VMA / HVA / 5-HIAA, lyophilised

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Intended use:

ClinTest® Standard Solutions are intended to check the daily performance of your HPLC - system. This Standard Solution is lyophilised and needs to be reconstituted before usage. Afterwards inject 20 µl directly into your analysing system. After separation, please compare the retention times of the peaks with those of earlier analyses. If you run the Standard Solution for the first time, you can compare your results with the chromatogram on the back of this data sheet. Besides retention times you are also able to check peak shape, integration and baseline with your standard chromatogram. For calibration, the ClinCal® Urine Calibrator, order no. 3013 has to be used.

Reconstitution:

Add exactly 1.0 ml of 0.2 M HCl to the vial and mix well.

Storage and stability:

Originally closed and stored at 2 - 8 °C, the lyophilised Standard Solution is stable for 36 months, but not longer than the expiration date printed on the label. After reconstitution the Standard Solution is stable for at least 48 hours when stored at 2 - 8 °C.

Pack size:

ClinTest® Standard Solution for VMA / HVA / 5-HIAA, lyophilised: 5 x 1 ml, **order no.: 3011**

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