ClinCal® - Calibrator

Serum Calibrator Set lyophilised / Serum-Kalibrator-Set lyophilisiert

FOR ANTIRETROVIRALS / FÜR ANTIRETROVIRALE MEDIKAMENTE

Intended use:

ClinCal® serum calibrators are used for calibration of the analysis-system. These lyophilised calibrators are based on human serum and are available as a calibrator set with 7 levels, including a blank. After reconstitution the calibrators have to be prepared like patient samples in one series of analyses.

Reconstitution:

Add exactly 0.5 ml of HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

Storage and stability:

Originally closed and stored at 2 - 8 °C, the lyophilised serum calibrators are stable for 36 months, but not longer than the expiration date printed on the label.

After reconstitution the stability of the analytes is:

- at least 5 hours when stored at 15 - 30 °C
- at least 1 day when stored at 2 - 8 °C
- at least 3 months when stored below -18 °C (avoid repeated freezing and thawing)

Notes:

The mean values of the analytes are chosen in a way that a linear calibration curve can be generated over a suitable range of concentration. According to quality assurance all ClinCal® calibrators have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume (CV) is < 1 %. The average residual moisture of this lot is 1.38 %.

Mean values:

The mean values have been established in independent reference laboratories with supervision of RECIPE, according to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK), with statistical methods.

Pack size:

ClinCal® Serum Calibrator Set
Level 0 - 6
7 x 1 x 0.5 ml, **order no.: MS4013**

Precautions:

The human serum which was used for manufacturing the calibrators was tested for the following infectious markers and found negative: HIV1/2-, HBV- and HCV-antibodies, Hepatitis B-surface antigen, HIV- and HCV-RNA, HBV-DNA (NAT). Nevertheless, the serum calibrators should be considered as potentially infectious and treated with appropriate care.

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REF MS4013
LOT 1326
2019-08

Level 0

This blank-level is manufactured from the same pool of human serum like the others, but it was not spiked with any antiretroviral drug. Therefore the concentration of each of the 16 analytes is zero. /

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Method of Analysis / Analysenmethode: LC-MS/MS

Please note: Dolutegravir and Ribavirine have been spiked additionally. Spike concentrations on request.

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RECIPE CHEMICALS + INSTRUMENTS GmbH