## ClinChek® - Control Plasma Control lyophilised / Kontrollplasma lyophilisiert

### FOR HOMOCYSTEINE / FÜR HOMOCYSTEIN

#### Intended use:

ClinChek® plasma controls are used for internal quality assurance in clinical-chemical laboratories. These lyophilised controls are based on human plasma and are available with mean values in the normal as well as in the pathological range. The controls have to be prepared like patient samples in one series of analyses.

#### Reconstitution:

Add exactly 3.0 ml of HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

If it is necessary to have a defined buffer-constitution of the control for further sample preparation (e.g. for pH-adjustment), you can also dissolve the lyophilised plasma control directly in the appropriate buffer.

#### Storage and stability:

Originally closed and stored at 2 - 8 °C, the lyophilised plasma control is stable for 36 months, but not longer than the expiration date printed on the label.

After reconstitution the stability of the analytes is:

- at least 12 hours when stored at 15 - 30 °C
- at least 7 days when stored at 2 - 8 °C
- at least 3 months when stored below -18 °C (avoid repeatedly freezing and thawing)

#### Notes:

Routine analyses of the controls are performed according to the testprocedure for patient samples. The concentration of the analyte is chosen in ranges where valid results can be obtained. According to quality assurance all ClinChek® controls have to pass strict quality-control procedures during manufacturing. REC-IPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume is max. +/- 0.25 %. The average residual moisture of this lot is 1.21 %.

#### Mean values:

The mean values and confidence intervals were established in additional independent reference laboratories with supervision of RECIFE, according to the guideline of the German Medical Council for quality assurance (RiliBäk).

#### Pack size:

ClinChek® Plasma Control  
Level I  
10 x 3 ml, **order no.: 23080**  
Level II  
10 x 3 ml, **order no.: 23081**  
Level I, II  
2 x 5 x 3 ml, **order no.: 23082**

#### Precautions:

The human plasma that was used for manufacturing the controls was tested for the following infectious markers and found negative: HIV 1/2- and HCV-antibodies, hepatitis B-surface antigen, HIV- and HCV-genome (PCR). Nevertheless, the plasma controls should be considered as potentially infectious and treated with appropriate care.

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RECIPE                      CHEMICALS + INSTRUMENTS GmbH