ClinChek[®] - Control

Serum Control lyophilised / Kontrollserum lyophilisiert

FOR ANTIPSYCHOTIC DRUGS TDM1 / FÜR PSYCHOPHARMAKA TDM1

Intended use:

ClinChek[®] serum controls are used for internal quality assurance in clinical-chemical laboratories. These lyophilised controls are based on human serum and are available with mean values in the therapeutic as well as in the elevated range. The controls have to be prepared like patient samples in one series of analyses.

Reconstitution:

Add exactly 5.0 ml of HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

If it is necessary to have a defined buffer-constitution of the control for further sample preparation (e.g. for pH-adjustment), you can also dissolve the lyophilised serum control directly in the appropriate buffer.

Storage and stability:

Originally closed and stored at 2 - 8 °C, the lyophilised serum controls are stable for 36 months, but not longer than the expiration date printed on the label.

After reconstitution the stability of the analytes is:

- at least 5 hours when stored at 15 - 30 °C
- at least 1 day when stored at 2 - 8 °C
- at least 3 months when stored below -18 °C (avoid repeatedly freezing and thawing)

Notes:

Routine analyses of the controls are performed according to the testprocedure for patient samples. The concentrations of the analytes are chosen in ranges where valid results can be obtained. According to quality assurance all ClinChek[®] controls have to pass strict quality-control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume is max. +/- 0.25 %. The average residual moisture of this lot is 0.91 %.

Mean values:

The mean values and confidence intervals have been established in independent reference laboratories with supervision of an institution that is recommended by the guideline of the German Medical Council for quality assurance (RiliBäk). For value assignment chromatographic methods were used.

Pack size:

ClinChek[®] Serum Control
Level I, II
2 x 5 x 5 ml, **order no.: 8702**

Precautions:

The human serum that was used for manufacturing the controls was tested for the following infectious markers and found negative: HIV 1/2-, HBc- and HCV-antibodies, hepatitis B-surface antigen, HIV-, HBV- and HCV-genome (PCR). Nevertheless, the serum controls should be considered as potentially infectious and treated with appropriate care.

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CHEMICALS + INSTRUMENTS GmbH Labortechnik