**ClinChek<sup>®</sup> - Control**  
**Serum Control lyophilised / Kontrollserum lyophilisiert**

<b>REF</b>	8850
<b>LOT</b>	919
	2012-05

**DRUGFREE****Intended use:**

ClinChek<sup>®</sup> serum controls are used for internal quality assurance in clinical-chemical laboratories. These lyophilised controls are based on calf serum. The controls have to be prepared like patient samples in one series of analyses.

**Reconstitution:**

Add exactly 10.0 ml of HPLCwater to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

If it is necessary to have a defined bufferconstitution of the control for further sample preparation (e.g. for pH-adjustment), you can also dissolve the lyophilised serum control directly in the appropriate buffer.

**Storage and stability:**

Originally closed and stored at 2 - 8 °C, the lyophilised serum control is stable for 36 months, but not longer than the expiration date printed on the label. After reconstitution the stability of the solution is:

- at least 5 hours  
when stored at 15 - 30 °C
- at least 24 hours  
when stored at 2 - 8 °C
- at least 10 days  
when stored below -18 °C  
(avoid repeatedly freezing and thawing)

**Notes:**

Routine analyses of the controls are performed according to the testprocedure for patient samples. According to quality assurance all ClinChek<sup>®</sup> controls have to pass strict quality-control procedures during manufacturing.

RECIPe guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume is max. +/- 0.25 %. The average residual moisture of this lot is 1.16%.

**Pack size:**

ClinChek<sup>®</sup> Serum Control  
Drugfree  
10 x 10 ml, **order no.: 8850**

**Precautions:**

The calf serum that was used for manufacturing the serum controls should be considered as potentially infectious and therefore the controls should be treated with appropriate care.

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