ClinChek[®] - Control

Urine Control lyophilised / Kontrollurin lyophilisiert

FOR 5-ALA, PBG / FÜR 5-ALA, PBG

Intended use:

ClinChek[®] urine controls are used for internal quality assurance in clinical-chemical laboratories. These lyophilised controls are based on human urine and are available with mean values in the normal as well as in the pathological range. The controls have to be prepared like patient samples in one series of analyses.

Reconstitution:

Add exactly 5.0 ml of HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use. If it is necessary to have a defined buffer-constitution of the control for further sample preparation (e.g. for pH-adjustment), you can also dissolve the lyophilised urine control directly in the appropriate buffer.

Storage and stability:

Originally closed and stored below -18 °C in the dark, the lyophilised urine control is stable for 36 months, but not longer than the expiration date printed on the label.

After reconstitution the stability of the analytes is when stored in the dark:

- at least 24 hours when stored at 15 - 30 °C
- at least 3 days when stored at 2 - 8 °C
- at least 3 months when stored below -18 °C

Notes:

Routine analyses of the controls are performed according to the testprocedure for patient samples. The concentrations of the analytes are chosen in ranges where valid results can be obtained. According to quality assurance all ClinChek[®] controls have to pass strict quality-control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume is max. +/- 0.25 %. The average residual moisture of this lot is 0.81 %.

Creatinine concentrations of these controls are: 0.37 g/l (level I) and 0.37 g/l (level II).

Mean values:

The mean values and confidence intervals were established in additional independent reference laboratories with supervision of RECIPE, according to the guideline of the German Medical Council (RililBäk).

Pack size:

ClinChek[®] Urine Control
Level I
10 x 5 ml, **order no.: 17280**
Level II
10 x 5 ml, **order no.: 17281**
Level I, II
2 x 5 x 5 ml, **order no.: 17282**

Precautions:

These controls were manufactured from human urine. Therefore the material has to be considered as potentially infectious and should be treated with appropriate care.

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CHEMICALS + INSTRUMENTS GmbH Labortechnik