ClinChek[®] - Control Urine Control lyophilised / Kontrollurin lyophilisiert

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Intended use:

ClinChek[®] urine controls are used for internal quality assurance in clinical-chemical laboratories. These lyophilised controls are based on human urine and are available with mean values in the normal as well as in the pathological range. The controls have to be prepared like patient samples in one series of analyses.

Reconstitution:

Add exactly 8.0 ml of 0.2 M HCl to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use. For analytics of 5-HIAA, we recommend to use 8.0 ml of 0.2 M acetic acid for reconstitution.

Storage and stability:

Originally closed and stored at 2 - 8 °C, the lyophilised urine control is stable for 36 months, but not longer than the expiration date printed on the label.

After reconstitution the stability of the analytes is:

- at least 12 hours
when stored at 15 - 30 °C
- at least 5 days
when stored at 2 - 8 °C
- at least 30 days
when stored below -18 °C
(avoid repeatedly freezing and thawing)

Notes:

Routine analyses of the controls are performed according to the testprocedure for patient samples. The concentrations of the analytes are chosen in ranges where valid results can be obtained. According to quality assurance all ClinChek[®] controls have to pass strict qualitycontrol procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume is max. +/- 0.25 %. The average residual moisture of this lot is 1.06 %.

Mean values:

The mean values and confidence intervals were established in additional independent reference laboratories with supervision of RECIPE, according to the guideline of the German Medical Council (RiliBäk) with statistic methods.

Pack size:

ClinChek[®] Urine Control
Level I
10 x 8 ml, **order no.: 8820**
ClinChek[®] Urine Control
Level II
10 x 8 ml, **order no.: 8821**
ClinChek[®] Urine Control
Level I, II
2 x 5 x 8 ml, **order no.: 8822**

Precautions:

These controls were manufactured from human urine. Therefore the material has to be considered as potentially infectious and should be treated with appropriate care.

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RECIPE CHEMICALS + INSTRUMENTS GmbH