ClinChek[®] - Control

Urine Control lyophilised / Kontrollurin lyophilisiert

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Intended use:

ClinChek[®] urine controls are used for internal quality assurance in the occupational and environmental toxicological laboratory as well as in the clinical-chemical one. These lyophilised controls are based on human urine and are available in two different ranges of concentration. After reconstitution, the controls have to be prepared like patient samples in one series of analyses.

Reconstitution:

Add exactly 10.0 ml of analytical grade or bidistilled water to the vial and close it well. Then dissolve the lyophilised material completely. When all material is dissolved, the solution is ready to use.

Storage and stability:

Originally closed and stored at 2 - 8 °C the lyophilised urine control is stable for 48 months, but not longer than the expiration date printed on the label.

After reconstitution the stability of the analytes is:

- at least 8 hours when stored at 15 - 30 °C
- at least 7 days when stored at 2 - 8 °C
- at least 1 month when stored below -18 °C (avoid repeatedly freezing and thawing)

If you intend to use reconstituted controls for the determination of mercury, zinc and aluminium, you will have a decrease of the value or contamination. Therefore we recommend to use freshly dissolved controls for these analytes.

Notes:

Routine analyses of the controls are performed according to the testprocedure for patient samples. The concentrations of the analytes are chosen in ranges where valid results can be obtained. According to quality assurance all ClinChek[®] controls have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume is max. +/- 0.25 %. The average residual moisture of this lot is 0.83 %.

Mean values:

The mean values and confidence intervals were established in independent reference laboratories with supervision of RECIPE, according to the guideline of the German Medical Council (RiliBäk) with statistic methods. For the determination of the trace elements, atomic absorption spectrometry (AAS) was used as well as ICP-MS. Even if you use different techniques of analysis your results should be within the given ranges of the data sheet.

Pack size:

ClinChek[®] Urine Control
Level I
10 x 10 ml, **order no.:** 8847

Level II
10 x 10 ml, **order no.:** 8848

Level I and II
2 x 5 x 10 ml, **order no.:** 8849

Precautions:

These controls were manufactured from human urine. Therefore the material has to be considered as potentially infectious and should be treated with appropriate care.

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CHEMICALS + INSTRUMENTS GmbH