## ClinChek® - Control

### Urine Control lyophilised / Kontrollurin lyophilisiert

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##### Intended use:

ClinChek® urine controls are used for internal quality assurance in clinical-chemical laboratories. These lyophilised controls are based on human urine and are available in three different levels of concentration. The controls have to be prepared like patient samples in one series of analyses.

##### Reconstitution:

Add exactly 3.0 ml of HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

##### Storage and stability:

Originally closed and stored at 2 - 8 °C, the lyophilised urine control is stable for 36 months, but not longer than the expiration date printed on the label. After reconstitution the stability of the analytes is:

- at least 12 hours when stored at 15 - 30 °C

- at least 7 days when stored at 2 - 8 °C
- at least 3 months when stored below -18 °C (avoid repeated freezing and thawing)

##### Notes:

Routine analyses of the controls are performed according to the test procedure for patient samples. The concentrations of the analytes are chosen in ranges where valid results can be obtained. According to quality assurance all ClinChek® controls have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume is max. +/- 0.25 %. The average residual moisture of this lot is 0.68 %.

Creatinine concentrations of the controls are: 0.22 g/l (level I), 0.51 g/l (level II) and 1.50 g/l (level III).

##### Mean values:

The mean values and confidence intervals were established in additional independent reference laboratories with supervision of RECIPE, according to the guideline of the German Medical Council (RiliBäk) with statistic methods.

##### Pack size:

ClinChek® Urine Control  
Level I  
5 x 3 ml, **order no.: MS3080**  
Level II  
5 x 3 ml, **order no.: MS3081**  
Level III  
5 x 3 ml, **order no.: MS3082**  
Level I, II, III  
3 x 2 x 3 ml, **order no.: MS3083**

##### Precautions:

These controls were manufactured from human urine. Therefore the material has to be considered as potentially infectious and should be treated with appropriate care.

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Technique of analysis: LC-MS/MS / Analysentechnik: LC-MS/MS

RECIPE CHEMICALS + INSTRUMENTS GmbH  
 Dessauerstr. 3 D-80992 Munich / Germany  
 Tel.: +49 / 89 / 54 70 81 - 0 Fax: +49 / 89 / 54 70 81 - 11  
 Internet: www.recipe.de e-mail: info@recipe.de