ClinChek® - Control

Hemolysate Control / Kontrollhämolysat

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Application:

ClinChek® Controls were used for internal quality control in clinical chemical laboratories. The Hemoglobin A1c-Control is a ready-to-use hemolysate of human blood with assigned values in the normal and elevated range.

Storage and stability:

Originally closed and stored at 2 - 8 °C, the material is stable for 24 months, but not longer than the expiration date printed on the label. After opening of the vials, the analytes are stable for 4 weeks, at 2 - 8 °C.

Notes:

Analysis of the controls does not require hemolysis, that means you can inject the material directly. During manufacturing, the ClinChek® Controls are subject to stringent quality control. RECIPE guarantees exactly the same composition and stability for each vial of a lot. Variation of filling volume is max. +/- 0.25 %.

Mean values:

The mean values and confidence intervals were established in additional independent reference laboratories with supervision of RECIPE, according to the guideline of the German Medical Council for quality assurance (RiliBäk). The Hb A1c values are calibrated in accordance with the DCCT (NGSP). Moreover, the controls are traceable to the IFCC-reference system.

Pack size:

ClinChek® Hemoglobin Control Level I
10 x 5 ml, **order no.: 11080**
Level II
10 x 5 ml, **order no.: 11081**
Level I, II
2 x 5 x 5 ml, **order no.: 11082**

Precautions:

The human whole blood that was used for manufacturing these controls was tested for the absence of hepatitis B-surface antigen and was HIV-negative. Nevertheless, the controls should be considered as potentially infectious and treated with appropriate care.

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$$\text{IFCC [mmol/mol]} = \text{NGSP [\%]} \times 10.93 - 23.5$$

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In Germany Hb A1c is subject to the guideline of the German Medical Council for quality assurance. For assessment of the internal quality control, please also refer to table B 1a of the guideline.

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CHEMICALS + INSTRUMENTS GmbH Labortechnik