# ClinChek® - Control

## Control-Hemolysate / Kontroll-Hämolysat

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#### Application:

ClinChek® Controls are used for internal quality control in clinical chemical laboratories. The Hemoglobin F/A2/S-Control is a ready-to-use hemolysate of human blood with assigned values in the elevated range.

#### Storage and stability:

Originally closed and stored at 2 - 8 °C, the material is stable for 36 months, but not longer than the expiration date printed on the label. After opening of the vials, the analytes are stable for 4 weeks, at 2 - 8 °C.

#### Notes:

Analysis of the control does not require hemolysis, that means you can inject the material directly in the analysing system. During manufacturing, the ClinChek® Controls are subject to stringent quality control. RECIPE guarantees exactly the same composition and stability for each vial of a lot. Variation of filling volume is max. +/- 0.25 %.

#### Mean values:

The mean values and confidence intervals were established in additional independent reference laboratories with supervision of RECIPE, according to the guideline of the German Medical Council for quality assurance (RiliBäk). The values for Hemoglobins F and A2 are traceable to the WHO-standards, NIBSC codes 85/616 and 89/666, respectively.

#### Pack size:

ClinChek® Control Hemolysate Level I  
5 x 3 ml, **order no.: 13183**

#### Precautions:

The human whole blood that was used for manufacturing these control was tested for the absence of hepatitis B-surface antigen and was HIV-negative. Nevertheless, the control should be considered as potentially infectious and treated with appropriate care.

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CHEMICALS + INSTRUMENTS GmbH      Labortechnik