## ClinCal® - Calibrator

**Serum Calibrator lyophilised, for Vitamin A and E /**

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The mean values have been assigned with the ClinRep® Complete Kits for Vitamin A and E (HPLC and UHPLC). The uncertainty of the mean values is  $\pm 1.7\%$  for Vitamin A and  $\pm 1.3\%$  for Vitamin E at a confidence level of 95%.  
**The mean values are traceable to NIST-SRM968e.**

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**Intended use:**

ClinCal® Serum Calibrators are used for calibration of the HPLC-system. After reconstitution this lyophilised calibrator has to be prepared like a patient sample.

**Reconstitution:**

Add exactly 3.0 ml HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

**Storage and stability:**

This product will be stable until the expiration date when stored unopened below  $-18\text{ °C}$  in the dark. After reconstitution the stability of the analytes is, when stored in the dark:

- 8 hours at  $15 - 30\text{ °C}$
- 3 days at  $2 - 8\text{ °C}$
- 90 days at below  $-18\text{ °C}$  (avoid repeated freezing and thawing)

The variation of the filling volume (CV) is  $< 1\%$ . The average residual moisture of this lot is  $0.84\%$ .

**Pack size:**

ClinCal® Serum Calibrator for Vitamin A and E: 1 x 3 ml, **order no.: 22013**

**Caution:**

The human serum which was used for manufacturing the calibrator was tested for the following infectious markers and found negative: HIV1/2- and HCV-antibodies, Hepatitis B-surface antigen, HIV- and HCV-RNA, HBV-DNA (NAT). Nevertheless, the serum calibrator should be considered as potentially infectious and treated with appropriate care.

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