**ClinCal® - Calibrator**  
**Serum Calibrator Set lyophilised / Serum-Kalibrator-Set lyophilisiert**

**FOR AMIODARONE / FÜR AMIODARON**

**Intended use:**

ClinCal® Serum Calibrators are used for calibration of the analysis-system. These lyophilised calibrators are based on calf serum and are available as a calibrator set with 4 levels, including a blank. After reconstitution the calibrators have to be prepared like patient samples in one series of analyses.

After reconstitution the stability of the analytes is:

- 8 hours when stored at 15 - 30 °C
- 5 days when stored at 2 - 8 °C
- 90 days when stored below -18 °C (avoid repeated freezing and thawing)

**Mean values:**

The mean values have been established in additional independent reference laboratories with supervision of RECIFE, according to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK), with statistical methods.

**Reconstitution:**

Add exactly 2.0 ml of HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

**Notes:**

The mean values of the analytes are chosen in a way that a linear calibration curve can be generated over a suitable range of concentration. According to quality assurance all ClinCal® Calibrators have to pass strict quality control procedures during manufacturing. RECIFE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume (CV) is < 1 %. The average residual moisture of this lot is 1.16 %.

**Pack size:**

ClinCal® Serum Calibrator Set Level 0 - 3  
 4 x 1 x 2 ml, **order no.: 9909**

**Storage and stability:**

This product will be stable until the expiration date when stored unopened below -18 °C.

**Precautions:**

This ClinCal® Calibrator Set was manufactured from pooled calfserum. It should be considered as potential infectious and treated with appropriate care.

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**REF** 9909  
**LOT** 545  
 2020-11

## Level 0

This blank-level was manufactured from the same pool of calf serum like the others, but it was not spiked with Amiodarone or Desethylamiodarone.

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RECIPE CHEMICALS + INSTRUMENTS GmbH