## ClinCal® - Calibrator

### Serum Calibrator Set lyophilised / Serum-Kalibrator-Set lyophilisiert

#### FOR STEROIDS / FÜR STEROIDE

##### Intended use:

ClinCal® Serum Calibrators are used for calibration of the analysis-system. These lyophilised calibrators are based on human serum and are available as a calibrator set with 7 levels, including a blank. After reconstitution the calibrators have to be prepared like patient samples in one series of analyses.

After reconstitution the stability of the analytes is:

- 10 days when stored at 2 - 8 °C
- 3 months when stored below -18 °C (avoid repeated freezing and thawing)

##### Notes:

The mean values of the analytes are chosen in a way that a linear calibration curve can be generated over a suitable range of concentration. According to quality assurance all ClinCal® Calibrators have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume (CV) is < 1 %. The average residual moisture of this lot is 1.04 %.

##### Reconstitution:

Add exactly 2.0 ml of HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

##### Storage and stability:

This product will be stable until the expiration date when stored unopened below -18 °C.

##### Mean values:

The mean values have been established in additional independent reference laboratories with supervision of RECIPE, according to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK), with statistical methods.

##### Pack size:

ClinCal® Serum Calibrator Set  
Level 0 - 6  
7 x 1 x 2 ml, **order no.: MS12013**

##### Precautions:

The human serum which was used for manufacturing the calibrators was tested for the following infectious markers and found negative: HIV1/2- and HCV-antibodies, Hepatitis B-surface antigen. Nevertheless, the serum calibrators should be considered as potentially infectious and treated with appropriate care.

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REF MS12013  
LOT 1226  
2019-06

## Level 0

The concentrations of the analytes are close to the lower limit of quantification of the analytical methods, which have been used for value assignment. Level 0 shall not be included in the calibration curve, rather shall be set as an „unknown“.

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## Method of Analysis / Analysenmethode: LC-MS/MS

### \* Please note / bitte beachten:

The mean values are traceable to **NIST-SRM971** and **NIST-SRM1950** /

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RECIPE CHEMICALS + INSTRUMENTS GmbH