ClinChek® - Control

Urine Control lyophilised / Kontrollurin lyophilisiert

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Intended use:

ClinChek® Urine Controls are used for internal quality assurance in the occupational and environmental toxicological laboratory as well as in the clinical-chemical one. These lyophilised controls are based on human urine and are available in two different ranges of concentration. After reconstitution the controls have to be prepared like patient samples in one series of analyses.

Reconstitution:

Add exactly 10.0 ml of analytical grade or double distilled water to the vial and close it well. Then dissolve the lyophilised material completely. When all material is dissolved, the solution is ready to use.

Storage and stability:

This product will be stable until the expiration date when stored unopened at 2 - 8 °C. After reconstitution the stability of the analytes is:

- 8 hours when stored at 15 - 30 °C
- 7 days when stored at 2 - 8 °C
- 30 days when stored below -18 °C (avoid repeated freezing and thawing)

If you intend to use reconstituted controls for the determination of mercury, zinc and aluminium, you will have a decrease of the value or contamination. Therefore we recommend to use freshly dissolved controls for these analytes.

Notes:

The concentrations of the analytes are chosen in ranges where valid results can be obtained. According to quality assurance all ClinChek® Controls have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume (CV) is < 1 %. The average residual moisture of this lot is 1.19 %.

Mean values:

The mean values and confidence intervals have been established in independent reference laboratories with supervision of RECIPE, according to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK), with statistical methods. For the determination of the trace elements ICP-MS was used. Even if you use different techniques of analysis, your results should be within the control ranges indicated on the data sheet.

Pack size:

ClinChek® Urine Control
Level I
10 x 10 ml, **order no.: 8847**
Level II
10 x 10 ml, **order no.: 8848**
Level I, II
2 x 5 x 10 ml, **order no.: 8849**

Precautions:

These controls were manufactured from human urine. Therefore the material has to be considered as potentially infectious and should be treated with appropriate care.

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*: For Gold we could not establish any mean values, due to the low number and the high imprecision of the available results (spike concentrations on request). /

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