## ClinChek® - Control Serum Control lyophilised / Kontrollserum lyophilisiert

### DRUGFREE

#### Intended use:

ClinChek® Serum Controls are used for internal quality assurance in clinical-chemical laboratories. These lyophilised controls are based on calf serum and are intended for use as a negative control (blank). After reconstitution the controls have to be prepared like patient samples in one series of analyses.

#### Reconstitution:

Add exactly 10.0 ml of HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

If it is necessary to have a defined bufferconstitution of the control for further sample preparation (e.g. for pH-adjustment), you can also dissolve the lyophilised serum control directly in the appropriate buffer.

#### Storage and stability:

This product will be stable until the expiration date when stored unopened at 2 - 8 °C. After reconstitution the stability of the serum is:

- 1 day when stored at 15 - 30 °C
- 5 days when stored at 2 - 8 °C
- 90 days when stored below -18 °C (avoid repeated freezing and thawing)

#### Notes:

According to quality assurance all ClinChek® Controls have to pass strict quality control procedures during manufacturing.

RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume (CV) is < 1 %. The average residual moisture of this lot is 2.76%.

#### Pack size:

ClinChek® Serum Control  
Drugfree  
10 x 10 ml, **order no.: 8850**

#### Precautions:

The calf serum which was used for manufacturing the serum controls was tested for the following viruses and found negative: BVD, IBR and PI3. Nevertheless, the serum controls should be considered as potentially infectious and treated with appropriate care.

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