ClinChek® - Control

Whole Blood Control lyophilised / Vollblutkontrolle lyophilisiert

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Intended use:

ClinChek® whole blood controls are used for internal quality assurance in clinical-chemical laboratories. These lyophilised controls are based on human whole blood and are available with mean values in the normal as well as in the elevated range.

Reconstitution:

Add exactly 0.5 ml of HPLC water to the vial. To dissolve the material completely, we recommend a so-called roller mixer. A treatment of 15 minutes on such an instrument is sufficient.

Procedure:

The reconstituted controls have to be pretreated the same as patient samples and analysed in accordance with instructions of the used analytic method.

Storage and stability:

Originally closed and stored at 2 - 8 °C, the lyophilised whole blood control is stable for 36 months, but not longer than the expiration date printed on the label.

After reconstitution the stability of the analyte is:

- at least 8 hours when stored at 15 - 30 °C
- at least 7 days when stored at 2 - 8 °C
- at least 30 days when stored below -18 °C (avoid repeated freezing and thawing)

Notes:

The concentrations of the analyte are chosen in ranges where valid results can be obtained. According to quality assurance all ClinChek® controls have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume is max. +/- 0.5 %. The average residual moisture of this lot is 0.89 %.

Total hemoglobin concentrations of the controls are: 14.5 g/dl (level I) and 14.3 g/dl (level II).

Mean values:

The mean values and confidence intervals have been established in independent reference laboratories with supervision of an institution that is recommended by the guideline of the German Medical Council for quality assurance (RiliBäk), with statistical methods. The Hb A1c values of the different methods are calibrated in accordance with the IFCC reference method. The metrological traceability is within the responsibility of the corresponding manufacturer of the Hb A1c-test.

Pack size:

ClinChek® Whole Blood Control
Level I, II
2 x 5 x 0.5 ml, **order no.: 11092**

Precautions:

The human whole blood which was used for manufacturing the controls was tested for the following infectious markers and found negative: HIV1, HIV2, HCV and hepatitis B-surface antigen. Nevertheless, the whole blood controls should be considered as potentially infectious and treated with appropriate care.

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RECIPE CHEMICALS + INSTRUMENTS GmbH