ClinChek® - Control **Hemolysate Control / Kontrollhämolysat** **FOR HEMOGLOBIN F/A2 / FÜR HÄMOGLOBIN F/A2**

REF	13180 - 13182
LOT	1516
	2022-12

Intended use:

ClinChek® Controls are used for internal quality control in clinical-chemical laboratories. The Hemoglobin F/A2-Control is a ready-to-use hemolysate of human blood with mean values in the normal and elevated range.

Storage and stability:

This product will be stable until the expiration date when stored unopened at 2 - 8 °C. After opening of the vials, the analytes are stable for 4 weeks, at 2 - 8 °C.

Notes:

Analysis of the controls does not require hemolysis, that means you can inject the material directly into your analysing system. According to quality assurance all ClinChek® Controls have to pass strict quality control procedures during manufacturing. RECIPE guarantees exactly the same composition and stability for each vial of a lot.

Mean values:

The mean values and confidence intervals have been established at RECIPE, according to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK), with statistical methods.

Pack size:

ClinChek® Hemolysate Control Level I, II
 2 x 5 x 3 ml, **order no.: 13182**

Precautions:

The human whole blood which was used for manufacturing these controls was tested for the absence of Hepatitis B-surface antigen and was HIV-negative. Nevertheless, the controls should be considered as potentially infectious and treated with appropriate care.

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