ClinTest® - Standard

Test Solution for CDT

for HPLC Assay in Serum

REF	21014
LOT	1269
	2021-06

Intended use:

ClinTest® Test Solutions are intended to check the daily performance of your HPLC - system. This Test Solution needs no preparation and can be injected directly into your analysing system (injection volume: 400 µl). In case of any turbidity, please centrifuge before injection. After separation, please compare the retention times of the peaks with those of earlier analyses. If you run the Test Solution for the first time, you can compare your results with the chromatogram below. Besides retention times you are also able to check peak shape, integration and baseline with your standard chromatogram.

Storage and stability:

The unopened Standard Solution, stored at 2 - 8 °C, can be used until the expiry date indicated on the product label. After use of the Standard Solution, the bottle must be closed tightly and stored at 2 - 8 °C, to ensure that it can be used until the expiry date indicated on the product label.


Pack size:

ClinTest® Test Solution for CDT: 1 x 15 ml, **order no.: 21014**

Precautions:

The human serum which was used for manufacturing the Test Solution was tested for the following infectious markers and found negative: HIV1/2- and HCV-antibodies as well as Hepatitis B-surface antigen. Nevertheless, the test solution should be considered as potentially infectious and treated with appropriate care.

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