
ClinChek® - Control

Hemoglobin Control lyophilised / Hämoglobinkontrolle lyophilisiert

8812-8814, 8816
2200
2023-05

Intended use:

ClinChek® Controls are used for internal quality assurance in clinical-chemical laboratories. These lyophilised controls are based on human erythrocytes and are available with mean values in three different ranges of concentration. After reconstitution the controls have to be prepared like patient samples in one series of analyses.

Reconstitution:

Add exactly 1.0 ml of reconstitution solution to the vial.

Caution:

The reconstitution solution, which is supplied with the controls, must be adapted to room temperature before use and can then be stored also at ambient temperature.

To dissolve the material completely, we recommend a so-called roller mixer. A treatment of 15 minutes on such an instrument is sufficient. In addition, we recommend the use of a vortex mixer.

Storage and stability:

This product will be stable until the expiration date when stored unopened at 2 - 8 °C. After reconstitution the stability of the analyte is:

- 1 day at 15 - 30 °C
- 2 days at 2 - 8 °C
- 90 days below -18 °C (avoid repeated freezing and thawing)

Notes:

The concentrations of the analyte are chosen in ranges where valid results can be obtained. According to quality assurance all ClinChek® Controls have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume (CV) is < 1 %.

Mean values:

The mean values and confidence intervals have been established at RECIPE, according to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK), with statistical methods.

Pack size:

ClinChek® Control
Level I (incl. reconstitution solution)
9 x 1 ml, **order no.: 8812**
Level II (incl. reconstitution solution)
9 x 1 ml, **order no.: 8813**
Level III (incl. reconstitution solution)
9 x 1 ml, **order no.: 8814**
Level I, II, III (incl. reconstitution solution)
3 x 3 x 1 ml, **order no.: 8816**

Precautions:

The human erythrocytes which were used for manufacturing the controls have been tested for the following infectious markers and found negative: HIV1/2-, HBV- and HCV-antibodies, Hepatitis B-surface antigen, HIV1- and HCV-RNA, HBV-DNA (NAT). Nevertheless, the controls should be considered as potentially infectious and treated with appropriate care.

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Technique of Analysis / Analysentechnik:

- (1) Cyanmethemoglobin method / Cyanmethämoglobin-Methode
- (2) 3-Wavelengths method according to Harboe / 3-Wellenlängen-Methode nach Harboe

IMPORTANT NOTE:

The Reconstitution Solution R (order no. 8815), which is included in the control packs, contains Triton X-100. As this substance is listed in Annex XIV (authorization list) of the REACH-Regulation (EC) No. 1907/2006, the use is prohibited as of 04. January 2021 with some exemptions.

With reference to articles 56 (3) and 56 (6), further use is allowed in the area of scientific research and development (SR&D), under controlled conditions and at concentrations below 0,1 % (w/w). This exemption for SR&D covers analytical applications in clinical laboratories, where Triton X-100 containing IVDs are used. Controlled conditions are considered to be given, provided that the use is limited to well trained laboratory professionals and that procedures for minimizing exposure / potential risks and professional waste disposal are established.

Further information is given on the Safety Data Sheet for the Reconstitution Solution R, which will be sent on demand. Please contact Customersupport@recipe.de.

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