ClinCal® - Calibrator

Serum Calibrator Set lyophilised / Serum-Kalibrator-Set lyophilisiert
NEUROLEPTICS / NEUROLEPTIKA

REF	MS9313
LOT	1368
	2020-11

Intended use:

ClinCal® Serum Calibrators are used for calibration of the analysis-system. These lyophilised calibrators are based on human serum and are available as a calibrator set with 4 levels, including a blank. After reconstitution the calibrators have to be prepared like patient samples in one series of analyses.

After reconstitution the stability of the analytes is:

- 1 day at 15 - 30 °C
- 7 days at 2 - 8 °C
- 90 days at below -18 °C (avoid repeated freezing and thawing)

Mean values:

The mean values have been established at RECIPE, according to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK), with statistical methods.

Pack size:

ClinCal® Serum Calibrator Set
 Level 0 - 3
 4 x 1 x 1 ml, **order no.: MS9313**

Reconstitution:

Add exactly 1.0 ml of HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

Notes:

The concentrations of the analytes are chosen in a way that a linear calibration curve can be generated over a suitable range of concentration. According to quality assurance all ClinCal® Calibrators have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume (CV) is < 1 %. The average residual moisture of this lot is 0.96 %.

Precautions:

The human serum which was used for manufacturing the calibrators was tested for the following infectious markers and found negative: HIV1/2- and HCV-antibodies, Hepatitis B-surface antigen, HIV- and HCV-RNA, HBV-DNA (NAT). Nevertheless, the serum calibrators should be considered as potentially infectious and treated with appropriate care.

Storage and stability:

This product will be stable until the expiration date when stored unopened at 2 - 8 °C.

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2020-11

Level 0

This blank-level was manufactured from the same pool of human serum like the others, but it was not spiked with any Neuroleptic Drug. /

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Method of Analysis / Analysenmethode: LC-MS/MS